

POWER FUNCTIONS



Power adjustable tilt



VR2 electronics standard with PS seating system



Expandable seat pan grows with you



Power seat elevator




Simple flip-up and removeable armrest

SHROUD COLORS



Crystal Blue




Vintage Wine




Metallic Green



Platinum Silver




Black Onyx



Chili Pepper Red



Popstar Pink



Titan Blue

SPECIFICATIONS

Maximum Speed	5 MPH	Power Seat Elevator	8"
Driving Range	15-19 Miles	Seat Widths	16" - 22" (Adjustable)
Max User Weight	300 lbs	Seat Depths	16" - 22" (Adjustable)
Ground Clearance	3.1"	Backrest Heights	22"
Battery Capacity	Group 34	Seat To Floor Height	16.75" - 19"
Width	24.25" - 28.5"	Electronics (standard)	VR2 Basic R-Net electronics optional
Length	43"	Alternative Drive Controls	Yes w/R-Net
Recommended Battery	Group 34	Effect of Driving	1/2 Speed @ 20° tilt angle
Gradient Capacity	10° Up 6° Dn		



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com